



## Decode a B Cell Fingerprint

A Viroclinics-DDL Publication - 2023



A Cerba Research Company

**WHITE**paper

# Contents


|  |           |
|--|-----------|
| Introduction                               | <b>3</b>  |
| How is B cell repertoire generated?        | <b>4</b>  |
| Viroclinics-DDL: One stop shop             | <b>5</b>  |
| BCR light/heavy chain analysis: An example | <b>6</b>  |
| Conclusion                                 | <b>12</b> |
| Appendix                                   | <b>13</b> |

# Introduction

The adaptive component of the immune system in humans is composed of a diverse population of B and T cells, which are paramount for the effective immune response upon exposure to a foreign antigens or mutated host ones (derived from pathogens, or mutated proteins expressed by cancer cells). In the case of B cells, antigens that are presented on the surface of an invading pathogen or mutated cell are recognized by a B cell receptor (BCR). Like a fingerprint being unique for any given person in the World, every BCR is unique for every B cell clone. Collectively, a pool of these clones generates a great variety of BCRs and antibodies, called the B-cell repertoire.

The application of next generation sequencing (NGS) technologies at present enables us to profile the B cell receptor repertoire with a high resolution. Profiling the BCR repertoire is typically done by assessing parameters such as clonality, diversity and richness of large number of B cells by sequencing their genetic code. Obtained genetic information is used to predict biophysical and other properties of a BCR or an antibody. This greatly enhances our ability to harness the potential of the selected clones for therapeutic, diagnostic or research applications.


B-cell repertoire profiling holds great promise to boost our efforts in fields of antibody development, immuno-therapeutics, vaccine development and disease monitoring. The ultra-deep high throughput sequencing we offer combined with state-of-the-art analyses developed by our partner is truly unparalleled.



# How is B cell repertoire generated?

B-cell receptors are essentially cell-surface-bound antibodies composed of a two light and a two heavy chains. Every B cell is expressing only one type of BCR and subsequently during the plasma cell stage, only one type of antibody. B cell receptors are generated through the process of genetic rearrangement of the large immunoglobulin gene segment. The diversity of BCRs arise from the process called somatic or V(D)J recombination. The gene loci of these receptors are organized as gene segment families comprising a variable (V) gene, a diversity (D) gene, a joining (J) gene and a constant (C) gene (**Figure 1**). V(D)J recombination presumes that rearrangement of these segments is followed by random nucleotide insertion and/or deletion at the junction points which creates even higher diversity of the B cell repertoire.


Both the heavy and light chain are composed of two domains: a variable (V) domain and a constant (C) domain. Each variable domain consists of three complementary-determining regions (CDRs). The CDR3 region is of specific interest as this region interacts with antigens and determines receptor specificity. The collection of all B-cells receptors genetically rearranged for different antigen specificity is known as the BCR repertoire.



**Figure 1:** Gene rearrangement in B cell receptor, heavy chain (light chain also depicted)

# Viroclinics-DDL: One stop shop

Viroclinics-DDL powers up with its partner ENPICOM to provide tailored, End to End solutions to its customers in decoding immune repertoires (**Figure 2**).



**Figure 2:** Schematic overview of the End-to-End workflow at Viroclinics-DDL


Viroclinics-DDL utilizes the Takara SMARTer BCR Profiling Kit that couples SMART (Switching Mechanism at 5' End of RNA Template) technology with an 5'RACE approach (**Appendix, Figure 1**), for unbiased amplification of BCR repertoire sequences. The repertoires are sequenced on Viroclinics-DDL's Illumina platforms and analyzed using the IGX platform developed by our partner ENPICOM. Viroclinics-DDL also support the development of tailor-made novel cellular and molecular assays, accompanied by advanced bioinformatic analysis.

# BCR light/heavy chain analysis: An example

## Quality control is essential to obtain useful data – RNA quality

Various factors may influence the quality of BCR profiling, for example the quality of samples being analyzed. The purpose of the example experiment outlined below is to determine what are the optimal sample amounts needed to obtain consistent, high-quality results.

Three separate isolations were performed on a human PBMC sample after which aliquots were made for each isolation to mimic biological and technical replicates. To address the question of how much input material is needed for successful BCR profiling, the aliquots were made with different inputs of RNA ranging from 50 to 800 ng.




To assess quantity and ensure quality of the isolated material, each sample is examined using the Bioanalyzer. RNA samples with a RNA integrity number (RIN) of 7 or higher are considered of sufficient quality to proceed with unbiased amplification of the BCR whole repertoire. All samples in the experiment met the quality criteria of RIN value of 7 or higher and were included in the BCR profiling assay (**Figure 3**).

**Figure 3:** All samples in the experiment had higher than required quality of RNA, reflected in the RIN value higher than 7.

## Quality control is essential to obtain useful data – Spike in controls to generate reproducible data

To ensure the quality of our workflow, control oligonucleotides (spike-in controls) are designed and spiked into the sample (**Appendix, Figure 10**). This provides a within-sample standard for between sample comparison and allows us to identify and correct for any amplification bias.

A spike-in control should never interfere with the sample itself because this will lead to a biased and unreliable result. Therefore, a minimal / maximal input concentration is determined. For each spike-in control concentration, the clone count fraction is determined (Figure 4). In the reactions spiked with 0.075 ng of spike-in control the fraction of reads used to sequence the spike-in control is 0.04, which is equal to 4% of the total reads. This amount is the upper limit in our runs.



**Figure 4:** Recovery of spike-in reads. Three concentrations of spike-in control are spiked in 200 ng PBMC RNA. Three samples with 0.075 ng spike-in control (Magenta), three samples with 0.15 ng spike-in control (green), three samples with 0.75 ng spike-in control (blue). Two reactions with maximum spike-in input, 7.125 ng, are the controls (purple). The three spike-in concentrations need respectively 4, 8, 15 and 100% of the total reads per sample.

## BCR Profiling: The analysis

The analysis of BCR repertoire sequencing data is not straightforward and requires different disciplines, including immunology and bioinformatics. The Viroclinics-DDL/ENPICOM team harbors the knowledge and expertise in all relevant areas to provide our customers with solid and understandable immune repertoire sequencing data. This can be either in a standard analysis format or tailored to address specific customer questions.

Our standard data analysis pipeline consists of the following elements:

- Raw/trimmed sequencing data quality
- Number of reads per UMI
- Total/unique clones
- Correlation analysis
- CDR3 length
- Richness
- Diversity
- V(D)J usage (top clones)
- AA sequence (top clones)
- Clonal distribution


Some of the standard reporting elements are described in more detail below using the human PBMC sample as an example.

## Number of total and unique clones recovered

For this experiment, an increasing amount of input RNA material was used to evaluate the range of nucleic acid input that provides accurate BCR profiling results. It is expected that with increasing nucleic acid input, an increasing number of BCR receptors (both total number of BCR receptors and unique BCR receptors) will be recovered. To evaluate reproducibility, technical replicates were included.

With an increasing concentration of RNA input we observe a steady increase in the number of total clones detected. An input of 400 ng resulted in the highest number of unique clones, although a significant part of unique receptors was already recovered using 50 ng of RNA input. The number of recovered total and unique clones is found to be highly similar between technical replicates.

**Conclusion:** As expected, the higher RNA input results in a higher number of receptors recovered. Using amounts between 200 ng and 400 ng of RNA input resulted in recovery of the highest number of unique receptors, therefore our recommendation is to use amounts in this range for successful BCR repertoire profiling. Finally, our results show that the workflow is highly reproducible in numbers of receptors recovered using technical replicates.



**Figure 5:** BCR receptors recovered using samples with increasing amount of RNA (range between 50 ng/μL and 800 ng/μL). The full bar represents the absolute number of BCR receptors recovered, the part in bold represent the unique BCR receptors recovered from the data.




## High level of similarity between technical replicates

To evaluate the reproducibility of our BCR profiling workflow in more depth, we evaluated similarity between technical replicates. The similarity in BCR clones identified between technical replicates is addressed by calculating the Ruzicka Similarity Coefficient. This coefficient takes clone count into consideration when assessing similarity between samples.

We observe high similarity between technical replicates, which is indicated in **Figure 6**. Similarity increases with an increasing amount of input material. As expected, 800 ng of input RNA represents the highest similarity (Upper right corner), whereas > 200 ng of input RNA acceptable similarity is observed.

**Conclusion:** Higher concentration of RNA input results in higher similarity between replicates. Based on similarity it is recommended to have at least 200 ng RNA as input.




**Figure 6:** Sample similarity plot: Technical replicates are plotted against each other and Ruzicka similarity coefficient was calculated. Similarity increases with an increased sample input amount (darker shades of green).

## Diversity & Richness

To assess diversity and richness as a function of RNA input amount, we plotted the number of unique clones recovered.

We observe that recovered diversity and richness steadily rose with an increase of the input RNA up to 200 ng. This was accompanied by high reproducibility between technical replicates.

**Conclusion:** The best recovery for both diversity and richness was obtained using an input amount of 200 ng of RNA. Considering other analyzed parameters, our recommendation is to use sample amount varying between 200 and 400 ng RNA. 800 ng of RNA input material results in a reduction of diversity and richness.




**Figure 7:** Diversity and richness as a function of the amount of the input RNA. Solid and reproducible recovery was observed using 200 and 400 ng of RNA.

## Top-10 B cell clones by CDR3 amino acid sequence

Analysis of the top clones was performed to examine the clonal distribution within the samples, depth of the reads per clone and the similarity in the number of reads between technical replicates in respect to the most abundant clones. The top-10 clones according to their CDR3 domain (amino acid sequence) are represented in **Figure 8**.

**Conclusion:** All five dilutions have a high overlap among the most abundant clones detected. For every top 10 clone, the number of reads recovered is highly similar between technical replicates (**Appendix, Table 1**). As expected, the number of reads per clone increases with the increasing amount of the input material.



**Figure 8:** Top 10 clones. Top clones are recovered from all five different input RNA amounts. Number of reads per clone was closely correlated with the RNA input amount.

# Conclusion

In this whitepaper we present the BCR sequencing workflow utilized at Viroclinics-DDL and discuss its features and benefits. Experimental setup of the presented case demonstrates our standard BCR repertoire sequencing workflow, the quality of generated data and the accuracy of our data analysis and interpretation.

## Take home messages:

- More RNA input results in a higher number of receptors recovered. 200 and 400 ng of RNA input yielded the highest recovery of unique receptors.
- Our results show that the BCR profiling workflow is highly reproducible in numbers of receptors recovered using technical replicates.
- Higher concentration of RNA input results in higher similarity between replicates. Based on similarity, it is recommended to have at least 200 ng of RNA as input. This also depends on the underlying scientific question and other parameters of interest.
- Diversity and richness are both highly similar between replicates. Best recovery is obtained with 200 ng RNA input. Our recommendation is to use between 200 to 400 ng of RNA for the optimal recovery of these parameters
- Finally, we examined clonal distribution for the most abundant clones according to amino acid sequence of the CDR3 domain. Same top clones are detected in all the samples, regardless of the amount of input material. Increasing amount of the RNA material results in the higher number of reads per clone.


In conclusion, our Viroclinics-DDL repertoire sequencing service provides reproducible high quality BCR profiling data. Optimal recovery of the unique clones, diversity and richness is obtained when using an RNA input ranging between 200 and 400 ng. Spike-in controls we use complete our end-to-end solution, a one-stop-shop for your immune profiling request.

If you want to know more about our **TCR** repertoire profiling services, check out the whitepaper: **Decode a T cell fingerprint** – available soon

# Appendix

## Technology in use


Takara Human BCR Profiling kit utilizes SMART technology (Switching Mechanism at 5' End of RNA Template) by combining NGS with a 5' RACE system for rapid amplification (**Figure 1**). In this approach, first strand cDNA synthesis starts from the 3' RNA poly A tail and it is dT-primed. During this step template switching reverse transcription reaction generates cDNAs with a sequence introduced by a template switching oligo attached to the 3' end of the cDNA (5' end of the transcript). Following, the first strand cDNA undergoes an enrichment process in which two rounds of gene-specific PCR amplifications are performed by using semi-nested PCR. The nested PCR reduces variability and allows for priming from the constant region of heavy or light chains. This method generates highly sensitive and reproducible B-cell repertoire profiling and allow to capture complete V(D)J variable regions of BCR transcripts.



**Figure 1:** Schematic overview of the Takara SMARTer BCR 5' RACE approach.

# Appendix

## Spike-in control design



**Figure 2:** Schematic overview of the IVTX1, V-gene/CDR3 (incl. AA sequence)/J-gene/C-gene.

| Sample ID | CDR3 AA Seq          | PBMC-50-1 | PBMC-50-2 | PBMC-100-1 | PBMC-100-2 | PBMC-200-1 | PBMC-200-2 | PBMC-400-1 | PBMC-400-2 | PBMC-800-1 | PBMC-800-2 |
|-----------|----------------------|-----------|-----------|------------|------------|------------|------------|------------|------------|------------|------------|
| 1         | CTRQGENGLAW          | 144       | 141       | 220        | 239        | 437        | 489        | 780        | 749        | 967        | 1048       |
| 2         | CAKGESYGDFMVDYYGMDVW | 92        | 116       | 216        | 180        | 390        | 379        | 651        | 599        | 753        | 889        |
| 3         | CARATTLGLAITYEYW     | 86        | 81        | 163        | 153        | 293        | 286        | 514        | 435        | 478        | 625        |
| 4         | CAKGTDVGGWDYW        | 80        | 92        | 175        | 155        | 293        | 306        | 531        | 463        | 602        | 728        |
| 5         | CSKHGNCRRGGDCYSDSW   | 80        | 94        | 148        | 176        | 313        | 340        | 594        | 571        | 730        | 848        |
| 6         | CVRDVGVMAMSHW        | 69        | 66        | 132        | 105        | 211        | 193        | 359        | 306        | 404        | 433        |
| 7         | CARLSDYGDYQRGFDPW    | 51        | 41        | 77         | NA         | 189        | 183        | 317        | 269        | 332        | 415        |
| 8         | CVRDTIIRGGYFDPW      | 51        | NA        | 78         | NA         | 160        | 172        | NA         | NA         | NA         | NA         |
| 9         | CAREIMATICFDYW       | 45        | 50        | 82         | 88         | NA         | 158        | 286        | 286        | 321        | 414        |
| 10        | CAREGHDLYHNLDW       | 44        | NA        | 79         | 89         | NA         | 168        | 277        | 247        | 307        | 357        |

Top 10 clones are stratified based on sample PBMC-50-1 \*

NA - Particular clone in respective dilution was not in its own top 10 \*

**Table 1:** Clonal distribution & the number of reads per clone.

Learn more about Viroclinics-DDL's  
Clinical Trial Operations Services on:  
[www.viroclinics.com/ctos](http://www.viroclinics.com/ctos)



# BCR RepSeq



## **Viroclinics-DDL**

Rotterdam Science Tower  
Marconistraat 16  
3029 AK Rotterdam  
The Netherlands

Tel. + 31 88 668 4787

[www.viroclinics.com](http://www.viroclinics.com)

[info@viroclinics.com](mailto:info@viroclinics.com)



A Cerba Research Company